

## Product Datasheet

### ACE-2 Recombinant Protein (orb1228246)

<b>Description</b>	ACE-2 Recombinant Protein
<b>Target</b>	Ace2
<b>Preservatives</b>	Supplied as a 0.2 um filtered solution of 25mM Tris-HCl, 150mM NaCl, 1mM ZnCl <sub>2</sub> , pH 7.5.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store below -20°C, stable for 6 months after receipt.
<b>Tag</b>	C-10 His tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant protein can be used for biological assays. For research use only.
<b>Protein Sequence</b>	Gln18-Thr740
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
<b>Source</b>	Mammalian Cells
<b>MW</b>	84.9KDa
<b>Uniprot ID</b>	<a href="#">Q8R0I0</a>
<b>NCBI</b>	<a href="#">Q8R0I0</a>
<b>Dilution Range</b>	This recombinant protein can be used for biological assays. For research use only.
<b>Expiration Date</b>	6 months from date of receipt.