



Product Datasheet

PRMT2 antibody (orb12283)


Description	Goat polyclonal to PRMT2.
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PRMT2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize all reported isoforms NP_001526.2;NP_001229793.1;NP_001229794.1;NP_001229795.1. Reported variants represent identical protein: NP_001526.2, NP_996845.1
Protein Sequence	DQPRRTTKYHSVILQ
Clonality	Polyclonal
MW	49.0; 37.2; 31.2; 31.2
Entrez	3275
Dilution Range	ELISA: 1:64000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate usi...



Western blot analysis of Mouse and Rat B...



Western blot analysis of Pig Brain lysat...