

---

## Product Datasheet

### Fc gamma RI Recombinant Protein (orb1228317)

**Description**

Fc gamma RI Recombinant Protein

**Target**

FCGR1A

**Preservatives**

Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4. It is not recommended to reconstitute to a concentration less than 100 ug/ml. Dissolve the lyophilized protein in ddH2O.

**Form/Appearance**

Lyophilized

**Storage**

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

**Tag**

C-6 His tag

**Note**

For research use only

**Application notes**

This recombinant protein can be used for biological assays. For research use only.

**Protein Sequence**

Gln16-Pro288

**Purity**

Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.

**Source**

Human Cells

**MW**

31.7 kD

**Uniprot ID**[P12314](#)**NCBI**[P12314](#)**Dilution Range**

This recombinant protein can be used for biological assays. For research use only.

**Expiration Date**

6 months from date of receipt.