

---

## Product Datasheet

### Peroxiredoxin-5 Recombinant Protein (orb1228383)

**Description**

Peroxiredoxin-5 Recombinant Protein

<b>Target</b>	PRDX5
<b>Preservatives</b>	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4. It is not recommended to reconstitute to a concentration less than 100 ug/ml. Dissolve the lyophilized protein in ddH2O.
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	N-6 His tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant protein can be used for biological assays. For research use only.
<b>Protein Sequence</b>	Met53-Leu214
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
<b>Source</b>	Human Cells
<b>MW</b>	17.9 kD
<b>Uniprot ID</b>	<a href="#">P30044</a>
<b>NCBI</b>	<a href="#">P30044</a>
<b>Dilution Range</b>	This recombinant protein can be used for biological assays. For research use only.
<b>Expiration Date</b>	6 months from date of receipt.