
Product Datasheet

Cardiotrophin-1 Recombinant Protein (orb1228720)

Description

Cardiotoxin-1 Recombinant Protein

Target	CTF1
Preservatives	Supplied as a 0.2 um filtered solution of 20mM PB \square 150mM NaCl,pH7.4. It is not recommended to reconstitute to a concentration less than 100 ug/ml.
Form/Appearance	Liquid
Storage	Store at -20°C, stable for 6 months after receipt. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Tag	N-6 His tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Met1-Ala201
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	E. coli
MW	22.6 kD
Uniprot ID	Q16619
NCBI	Q16619
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.