
Product Datasheet

Resistin Recombinant Protein (orb1228762)

Description

Resistin Recombinant Protein

Target	RETN
Preservatives	Lyophilized from a 0.2 um filtered solution of 20mM Acetic acid,pH 3.0. It is not recommended to reconstitute to a concentration less than 100 ug/ml. Dissolve the lyophilized protein in ddH2O.
Form/Appearance	Lyophilized
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6 His tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Lys19-Pro108
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	E. coli
MW	10.6 kD
Uniprot ID	Q9HD89
NCBI	Q9HD89
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.