

Product Datasheet

UBE2Z Recombinant Protein (orb1228777)

Description	UBE2Z Recombinant Protein
Target	UBE2Z
Preservatives	Supplied as a 0.2 um filtered solution of 50mM HEPES, 50mM NaCl, pH 8.0. It is not recommended to reconstitute to a concentration less than 100 ug/ml.
Form/Appearance	Liquid
Storage	Store at -20°C, stable for 6 months after receipt. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Tag	N-6 His tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Met1-Val246
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	E. coli
MW	30.2 kD
Uniprot ID	Q9H832
NCBI	Q9H832
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.