
Product Datasheet

Autophagy protein 10 Recombinant Protein (orb1228928)

Description

Autophagy protein 10 Recombinant

Target	ATG10
Preservatives	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 8.0. It is not recommended to reconstitute to a concentration less than 100 ug/ml.
Form/Appearance	Liquid
Storage	Store at -20°C, stable for 6 months after receipt. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Tag	C-6 His, N-T7 tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Met1-Thr190
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	E. coli
MW	24.5 kD
Uniprot ID	Q9H0Y0
NCBI	Q9H0Y0
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.