

---

## Product Datasheet

### Autophagy protein 3 Recombinant Protein (orb1228930)

**Description**

Autophagy protein 3 Recombinant

<b>Target</b>	ATG3
<b>Preservatives</b>	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.2. It is not recommended to reconstitute to a concentration less than 100 ug/ml.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store at -20°C, stable for 6 months after receipt. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
<b>Tag</b>	C-6 His, N-T7 tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant protein can be used for biological assays. For research use only.
<b>Protein Sequence</b>	Met1-His287
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
<b>Source</b>	E. coli
<b>MW</b>	35 kD
<b>Uniprot ID</b>	<a href="#">Q9NT62</a>
<b>NCBI</b>	<a href="#">Q9NT62</a>
<b>Dilution Range</b>	This recombinant protein can be used for biological assays. For research use only.
<b>Expiration Date</b>	6 months from date of receipt.