

---

## Product Datasheet

### CD300a Recombinant Protein (orb1229007)

**Description**

CD300a Recombinant Protein

<b>Target</b>	CD300A
<b>Preservatives</b>	Lyophilized from a 0.2 um filtered solution of 20mM Tris, 150mM NaCl, 5% Trehalose, pH 8.0. It is not recommended to reconstitute to a concentration less than 100 ug/ml. Dissolve the lyophilized protein in ddH <sub>2</sub> O.
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	C-Fc-6 His tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant protein can be used for biological assays. For research use only.
<b>Protein Sequence</b>	Leu18-Gln178
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
<b>Source</b>	Human Cells
<b>MW</b>	45.4 kD
<b>Uniprot ID</b>	<a href="#">Q9UGN4</a>
<b>NCBI</b>	<a href="#">Q9UGN4</a>
<b>Dilution Range</b>	This recombinant protein can be used for biological assays. For research use only.
<b>Expiration Date</b>	6 months from date of receipt.