

## Product Datasheet

### UBE2B Recombinant Protein (orb1229238)

<b>Description</b>	UBE2B Recombinant Protein
<b>Target</b>	UBE2B
<b>Preservatives</b>	Supplied as a 0.2 um filtered solution of 20mM PB,150mM NaCl,pH7.2. It is not recommended to reconstitute to a concentration less than 100 ug/ml.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store at -20°C, stable for 6 months after receipt. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
<b>Tag</b>	C-6 His tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant protein can be used for biological assays. For research use only.
<b>Protein Sequence</b>	Met1-Ser152
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
<b>Source</b>	Human Cells
<b>MW</b>	18.3 kD
<b>Uniprot ID</b>	<a href="#">P63146</a>
<b>NCBI</b>	<a href="#">P63146</a>
<b>Dilution Range</b>	This recombinant protein can be used for biological assays. For research use only.
<b>Expiration Date</b>	6 months from date of receipt.