

Product Datasheet

ACE-2 Recombinant Protein (orb1229507)

Description	ACE-2 Recombinant Protein
Target	ACE2
Preservatives	Supplied as a 0.2 um filtered solution of 20mMTris,300mMNaCl,1mMZnCl2,10%Glycerol,PH7.4 .
Form/Appearance	Liquid
Storage	Store at -20°C, stable for 6 months after receipt. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Tag	C-6 His tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Gln18-Ser740
Purity	Greater than 95% as determined by reducing SDS- PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	Human Cells
MW	84.63 kD
Uniprot ID	Q9BYF1
NCBI	Q9BYF1
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.