

Product Datasheet

MMP-3 Recombinant Protein (orb1229536)

Description	MMP-3 Recombinant Protein
Target	MMP3
Preservatives	Supplied as a 0.2 um filtered solution of 20mM TrisHCl, 150mM NaCl, 0.05% Brij35, 10% Glycerol, pH 7.5. It is not recommended to reconstitute to a concentration less than 100 ug/ml.
Form/Appearance	Liquid
Storage	Store at -20°C, stable for 6 months after receipt. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Tag	C-6 His tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Tyr18-Cys477
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	Human Cells
MW	53.26 kD
Uniprot ID	P08254
NCBI	P08254
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.