

Product Datasheet

OTUB2 Recombinant Protein (orb1229655)

Description	OTUB2 Recombinant Protein
Target	OTUB2
Preservatives	Supplied as a 0.2 um filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, pH 7.5 . It is not recommended to reconstitute to a concentration less than 100 ug/ml.
Form/Appearance	Liquid
Storage	Store at -20°C, stable for 6 months after receipt. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Tag	N-GST tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Met1-His234
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	E. coli
MW	53.6 kD
Uniprot ID	Q96DC9
NCBI	Q96DC9
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.