
Product Datasheet

IL-1F10 Recombinant Protein (orb1229803)

Description

IL-1F10 Recombinant Protein

Target	IL1F10
Preservatives	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 6.0. It is not recommended to reconstitute to a concentration less than 100 ug/ml. Dissolve the lyophilized protein in ddH2O.
Form/Appearance	Lyophilized
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	Tag Free
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Met1-Trp152
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	E. coli
MW	16.9 kD
Uniprot ID	Q8WWZ1
NCBI	Q8WWZ1
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.