
Product Datasheet

IFN gamma Recombinant Protein (orb1229849)

Description

IFN gamma Recombinant Protein

Target	IFNG
Preservatives	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.0. It is not recommended to reconstitute to a concentration less than 100 ug/ml. Dissolve the lyophilized protein in ddH2O.
Form/Appearance	Lyophilized
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	Tag Free
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Gln24-Gln166
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	E. coli
MW	16.88 kD
Uniprot ID	P01579
NCBI	P01579
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.