
Product Datasheet

IL-11 Recombinant Protein (orb1229856)

Description

IL-11 Recombinant Protein

| | |
|--------------------------|---|
| Target | IL11 |
| Preservatives | Lyophilized from a 0.2 um filtered solution of 20mM PB, 2% Glycine, pH 7.2. It is not recommended to reconstitute to a concentration less than 100 ug/ml. Dissolve the lyophilized protein in ddH2O. |
| Form/Appearance | Lyophilized |
| Storage | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. |
| Tag | Tag Free |
| Note | For research use only |
| Application notes | This recombinant protein can be used for biological assays. For research use only. |
| Protein Sequence | Gly23-Leu199 |
| Purity | Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test. |
| Source | P. pastoris |
| MW | 19 kD |
| Uniprot ID | P20809 |
| NCBI | P20809 |
| Dilution Range | This recombinant protein can be used for biological assays. For research use only. |
| Expiration Date | 6 months from date of receipt. |