

Product Datasheet

IL-35 Recombinant Protein (orb1229883)

Description	IL-35 Recombinant Protein
Target	EBI3
Preservatives	Lyophilized from 0.2um-filtered solution in PBS. Reconstitute at 100 µg/ml in sterile PBS.
Form/Appearance	Lyophilized
Storage	Stable for at least 1 year after receipt when stored at -20°C. Working aliquots are stable for up to 3 months when stored at -20°C.
Tag	Fc Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	The human IL-35 complex composed of the Ebi3 subunit (aa 21-229) and the IL-12a (p35) subunit (aa 57-253) is fused through a polypeptide linker to the Fc portion of human IgG1.
Purity	> 98% (SDS-PAGE). Endotoxin level is less than 0.06EU/ µg protein (LAL test; Lonza).
Source	CHO cells
Uniprot ID	Q4QNY4
NCBI	NP_005746 , NP_000873
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.