

Product Datasheet

IL-22 Recombinant Protein (orb1229898)

| | |
|--------------------------|---|
| Description | IL-22 Recombinant Protein |
| Target | IL22 |
| Preservatives | Lyophilized from 0.2um-filtered solution in PBS. Reconstitute at 100 µg/ml in sterile PBS. |
| Form/Appearance | Lyophilized |
| Storage | Stable for at least 1 year after receipt when stored at -20°C. Working aliquots are stable for up to 3 months when stored at -20°C. |
| Tag | Fc Tag |
| Note | For research use only |
| Application notes | This recombinant proteins is for research use only. |
| Protein Sequence | The extracellular domain of mouse IL-22 (aa 34-179) is fused to the N-terminus of the Fc region of mouse IgG2a. |
| Purity | > 98% (SDS-PAGE). Endotoxin level is less than 0.06EU/ µg protein (LAL test; Lonza). |
| Source | CHO cells |
| Uniprot ID | P63159 |
| NCBI | NP_058667 |
| Dilution Range | This recombinant proteins is for research use only. |
| Expiration Date | 6 months from date of receipt. |