

Product Datasheet

IL-21R non-lytic Recombinant Protein (orb1229905)

Description	IL-21R non-lytic Recombinant Protein
Target	Il21r
Preservatives	Lyophilized from 0.2um-filtered solution in PBS. Reconstitute at 100 µg/ml in sterile PBS.
Form/Appearance	Lyophilized
Storage	Stable for at least 1 year after receipt when stored at -20°C. Working aliquots are stable for up to 3 months when stored at -20°C.
Tag	Fc Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	The extracellular domain of mouse IL-21R (aa 20-234) is fused to the N-terminus of the Fc region of a mutant mouse IgG2a.
Purity	> 98% (SDS-PAGE). Endotoxin level is less than 0.06EU/ µg protein (LAL test; Lonza).
Source	CHO cells
Uniprot ID	Q9EPC2
NCBI	AAL82636
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.