

Product Datasheet

IL-21 mutant Recombinant Protein (orb1229907)

Description	IL-21 mutant Recombinant Protein
Target	IL21
Preservatives	Lyophilized from 0.2um-filtered solution in PBS. Reconstitute at 100 µg/ml in sterile PBS.
Form/Appearance	Lyophilized
Storage	Stable for at least 1 year after receipt when stored at -20°C. Working aliquots are stable for up to 3 months when stored at -20°C.
Tag	His Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	Human IL-21 (aa 32-162) with mutations at Gln145Asp and Gln150Asp is fused at the N-terminus to a His-tag.
Purity	> 90% (SDS-PAGE). Endotoxin level is less than 1EU/ µg protein (LAL test; Lonza).
Source	E. coli
Uniprot ID	Q9UNE0
NCBI	NP_068575
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.