

Product Datasheet

IL-15 mutant Recombinant Protein (biotin) (orb1229913)

Description	IL-15 mutant Recombinant Protein (biotin)
Target	IL15
Preservatives	In PBS and 0.09% sodium azide.
Form/Appearance	Lyophilized
Storage	Stable for at least 6 months after receipt when stored at +4°C.
Tag	Fc Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	The extracellular domain of human IL-15 (aa 49-162) including mutations at Q149D and Q156D is fused to the N-terminus of the Fc region of mouse IgG2a.
Purity	> 98% (SDS-PAGE). Endotoxin level is less than 0.06EU/ µg protein (LAL test; Lonza).
Source	CHO cells
Uniprot ID	O88516
NCBI	NP_751914
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.