

Product Datasheet

IL-10 Recombinant Protein (orb1229920)

Description	IL-10 Recombinant Protein
Target	IL10
Preservatives	Lyophilized from 0.2um-filtered solution in PBS (pH 7.4). Reconstituted in sterile H2O not less than 100 µg/ml, which can then be further diluted in other aqueous solutions.
Form/Appearance	Lyophilized
Storage	Stable for at least 1 year after receipt when stored at -20°C. Working aliquots are stable for up to 3 months when stored at -20°C.
Tag	His Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	Human IL-10 (aa 19-178) is fused at the N-terminus to a His-tag.
Purity	> 95% (SDS-PAGE). Endotoxin level is less than 0.1EU/ µg protein (LAL test; Lonza).
Source	E. coli
Uniprot ID	Q8TDQ0
NCBI	NP_000563
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.