

## Product Datasheet

### IL-2 non-lytic Recombinant Protein (orb1229948)

<b>Description</b>	IL-2 non-lytic Recombinant Protein
<b>Target</b>	IL2
<b>Preservatives</b>	Lyophilized from 0.2um-filtered solution in PBS. Reconstitute at 100 µg/ml in sterile PBS.
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Stable for at least 1 year after receipt when stored at -20°C. Working aliquots are stable for up to 3 months when stored at -20°C.
<b>Tag</b>	Fc Tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant proteins is for research use only.
<b>Protein Sequence</b>	The extracellular domain of mouse IL-2 (aa 21-169) is fused to the N-terminus of the Fc region of a mutant mouse IgG2a.
<b>Purity</b>	> 98% (SDS-PAGE). Endotoxin level is less than 0.06EU/ µg protein (LAL test; Lonza).
<b>Source</b>	NS1 cells
<b>Uniprot ID</b>	<a href="#">Q9GZX6</a>
<b>NCBI</b>	<a href="#">NP_032392</a>
<b>Dilution Range</b>	This recombinant proteins is for research use only.
<b>Expiration Date</b>	6 months from date of receipt.