

Product Datasheet

CX3CL1 Recombinant Protein (orb1229963)

Description	CX3CL1 Recombinant Protein
Target	CX3CL1
Preservatives	Lyophilized from 0.2um-filtered solution in PBS. Reconstitute at 100 µg/ml in sterile PBS.
Form/Appearance	Lyophilized
Storage	Stable for at least 1 year after receipt when stored at -20°C. Working aliquots are stable for up to 3 months when stored at -20°C.
Tag	Fc Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	The extracellular domain of human CX3CL1 (aa 25-100) is fused to the N-terminus of the Fc region of human IgG1.
Purity	> 98% (SDS-PAGE). Endotoxin level is less than 0.06EU/ µg protein (LAL test; Lonza).
Source	CHO cells
Uniprot ID	P18893
NCBI	AAH01163
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.