

Product Datasheet

CD137 Recombinant Protein (orb1230006)

Description	CD137 Recombinant Protein
Target	TNFRSF9
Preservatives	Lyophilized from a concentrated sterile solution containing 50mM Tris-HCl buffer (pH 8.0) and 500mM NaCl. Reconstitute in sterile water not less than 100 µg/ml, which can then be further diluted to other aqueous solutions.
Form/Appearance	Lyophilized
Storage	Stable for at least 1 year after receipt when stored at -20°C. Working aliquots are stable for up to 3 months when stored at -20°C.
Tag	His Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	Human CD137 (aa 19-184) is fused at the C-terminus to a His-tag.
Purity	≥95% (SDS-PAGE). Endotoxin level is less than 0.1EU/ µg protein (LAL test; Lonza).
Source	E. coli
Uniprot ID	P01588
NCBI	AAH06196
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.