

Product Datasheet

CD80 Recombinant Protein (orb1230016)

Description	CD80 Recombinant Protein
Target	CD80
Preservatives	Lyophilized from 0.2um-filtered solution in PBS. Reconstitute with 100 µl sterile water. Add 1X PBS to the desired protein concentration.
Form/Appearance	Lyophilized
Storage	Stable for at least 1 year after receipt when stored at -20°C. Working aliquots are stable for up to 3 months when stored at -20°C.
Tag	Fc Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	The extracellular domain of human CD80 (aa 35-242) is fused to the N-terminus of the Fc region of mouse IgG2a.
Purity	> 98% (SDS-PAGE). Endotoxin level is less than 0.06EU/ µg protein (LAL test; Lonza).
Source	CHO cells
Uniprot ID	Q9WVS0
NCBI	NP_005182
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.