

## Product Datasheet

# NMNAT3 highly active Recombinant Protein (orb1230060)

<b>Description</b>	NMNAT3 highly active Recombinant Protein
<b>Target</b>	NMNAT3
<b>Preservatives</b>	In 50mM TRIS-HCl, pH 8.0, containing 300mM sodium chloride, 1mM DTT and 10% glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Working aliquots are stable for up to 3 months when stored at -20°C.
<b>Tag</b>	His Tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant proteins is for research use only.
<b>Protein Sequence</b>	Human full-length NMNAT3 (aa 1-252) is fused at the N-terminus to a His-tag.
<b>Purity</b>	> 95% (SDS-PAGE)
<b>Source</b>	E. coli
<b>Uniprot ID</b>	<a href="#">P333681</a>
<b>NCBI</b>	<a href="#">Q96T66</a>
<b>Dilution Range</b>	This recombinant proteins is for research use only.
<b>Expiration Date</b>	6 months from date of receipt.