

## Product Datasheet

# NF-kB p50 highly active Recombinant Protein (orb1230066)

<b>Description</b>	NF-kB p50 highly active Recombinant Protein
<b>Target</b>	NFKB1
<b>Preservatives</b>	In 50mM TRIS-HCl, pH 7.5, containing 100mM sodium chloride, 0.2% NP-40, 50-100mM imidazole and 10% glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	Lot dependent (approx. 0.2mg/ml)
<b>Storage</b>	Stable for at least 6 months after receipt when stored at -80°C.
<b>Tag</b>	His Tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant proteins is for research use only.
<b>Protein Sequence</b>	Human NF-kappaB (p50) is fused to a His-tag.
<b>Purity</b>	> 95% (SDS-PAGE)
<b>Source</b>	Sf21 cells
<b>Uniprot ID</b>	<a href="#">P26842</a>
<b>NCBI</b>	<a href="#">P19838</a>
<b>Dilution Range</b>	This recombinant proteins is for research use only.
<b>Expiration Date</b>	6 months from date of receipt.