

Product Datasheet

NF-kB p65 highly active Recombinant Protein (orb1230067)

Description	NF-kB p65 highly active Recombinant Protein
Target	RELA
Preservatives	In 50mM TRIS-HCl, pH 7.5, containing 100mM sodium chloride, 0.2% NP-40, 50mM imidazole and 10% glycerol.
Form/Appearance	Liquid
Concentration	Lot dependent (0.2-1.5mg/ml)
Storage	Stable for at least 6 months after receipt when stored at -80°C.
Tag	His Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	Human NF-kappaB (p65) is fused to a His-tag.
Purity	> 95% (SDS-PAGE)
Source	Sf21 cells
Uniprot ID	P26842
NCBI	Q04206
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.