

Product Datasheet

IL-38 Recombinant Protein (orb1230126)

Description	IL-38 Recombinant Protein
Target	Il1f10
Preservatives	Contains PBS + 0.5% Trehalose. Reconstitute with 100 µl sterile water.
Form/Appearance	Lyophilized
Concentration	0.1mg/ml after reconstitution.
Storage	Stable for at least 6 months after receipt when stored at -20°C.
Tag	DDDDK Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	Mouse IL-38 [IL-1F10] (aa 3-152) is fused at the C-terminus to a DDDDK-tag.
Purity	> 95% (SDS-PAGE). Endotoxin level is less than 0.01EU/ µg purified protein (LAL test; Lonza).
Source	E. coli
Uniprot ID	Q9Y275
NCBI	Q8R459
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.