

## Product Datasheet

### IL-36 gamma Recombinant Protein (orb1230127)

<b>Description</b>	IL-36 gamma Recombinant Protein
<b>Target</b>	Il1f9
<b>Preservatives</b>	Contains PBS + 0.5% Trehalose. Reconstitute with 100 µl sterile water.
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	0.1mg/ml after reconstitution.
<b>Storage</b>	Stable for at least 6 months after receipt when stored at -20°C.
<b>Tag</b>	Tag Free
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant proteins is for research use only.
<b>Protein Sequence</b>	Mouse IL-36gamma (aa Gly13-Ser164) is untagged.
<b>Purity</b>	> 95% (SDS-PAGE). Endotoxin level is less than 0.01EU/ µg purified protein (LAL test; Lonza).
<b>Source</b>	E. coli
<b>Uniprot ID</b>	<a href="#">P49771</a>
<b>NCBI</b>	<a href="#">Q8R460</a>
<b>Dilution Range</b>	This recombinant proteins is for research use only.
<b>Expiration Date</b>	6 months from date of receipt.