

Product Datasheet

IL-36 beta Recombinant Protein (orb1230128)

Description	IL-36 beta Recombinant Protein
Target	Il1f8
Preservatives	Contains PBS + 0.5% Trehalose. Reconstitute with 100 µl sterile water.
Form/Appearance	Lyophilized
Concentration	0.1mg/ml after reconstitution.
Storage	Stable for at least 6 months after receipt when stored at -20°C.
Tag	Tag Free
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	Mouse IL-36beta (aa Ser31-Lys183) is untagged.
Purity	> 95% (SDS-PAGE). Endotoxin level is less than 0.01EU/ µg purified protein (LAL test; Lonza).
Source	E. coli
Uniprot ID	Q9NZH7
NCBI	Q9D6Z6
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.