

Product Datasheet

IL-36 alpha Recombinant Protein (orb1230129)

Description	IL-36 alpha Recombinant Protein
Target	Il1f6
Preservatives	Contains PBS. Reconstitute with 100 µl sterile water.
Form/Appearance	Lyophilized
Concentration	0.1mg/ml after reconstitution.
Storage	Stable for at least 6 months after receipt when stored at -20°C.
Tag	Tag Free
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	Mouse IL-36alpha (aa Arg8-His160) is untagged.
Purity	> 95% (SDS-PAGE). Endotoxin level is less than 0.01EU/ µg purified protein (LAL test; Lonza).
Source	E. coli
Uniprot ID	Q8NAU1
NCBI	Q9JLA2
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.