

## Product Datasheet

### BCMA Recombinant Protein (orb1230144)

<b>Description</b>	BCMA Recombinant Protein
<b>Target</b>	TNFRSF17
<b>Preservatives</b>	Contains PBS. Reconstitute with 50 $\mu$ l sterile water.
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	1mg/ml after reconstitution.
<b>Storage</b>	Stable for at least 6 months after receipt when stored at -20°C.
<b>Tag</b>	Fc Tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant proteins is for research use only.
<b>Protein Sequence</b>	The extracellular domain of human BCMA (aa 2-54) is fused at the C-terminus to the Fc portion of human IgG1.
<b>Purity</b>	> 95% (SDS-PAGE). Endotoxin level is less than 0.01EU/ $\mu$ g purified protein (LAL test; Lonza).
<b>Source</b>	HEK293 cells
<b>Uniprot ID</b>	<a href="#">P15018</a>
<b>NCBI</b>	<a href="#">Q02223</a>
<b>Dilution Range</b>	This recombinant proteins is for research use only.
<b>Expiration Date</b>	6 months from date of receipt.