

Product Datasheet

TACI Recombinant Protein (orb1230145)

Description	TACI Recombinant Protein
Target	TNFRSF13B
Preservatives	Contains PBS. Reconstitute with 50 µl sterile water.
Form/Appearance	Lyophilized
Concentration	1mg/ml after reconstitution.
Storage	Stable for at least 6 months after receipt when stored at -20°C.
Tag	Fc Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	The extracellular domain of human TACI (aa 2-118) is fused at the C-terminus to the Fc portion of human IgG1.
Purity	> 90% (SDS-PAGE). Endotoxin level is less than 0.01EU/ µg purified protein (LAL test; Lonza).
Source	HEK293 cells
Uniprot ID	O35305
NCBI	O14836
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.