

Product Datasheet

TNF-R1 Recombinant Protein (orb1230150)

Description	TNF-R1 Recombinant Protein
Target	TNFRSF1A
Preservatives	Contains PBS. Reconstitute with 50 μ l sterile water.
Form/Appearance	Lyophilized
Concentration	1mg/ml after reconstitution.
Storage	Stable for at least 6 months after receipt when stored at -20°C.
Tag	Fc Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	The extracellular domain of human TNF-R1 (aa 22-211) is fused to the Fc portion of human IgG1.
Purity	> 95% (SDS-PAGE). Endotoxin level is less than 0.01EU/ μ g purified protein (LAL test; Lonza).
Source	HEK293 cells
Uniprot ID	P25942
NCBI	P19438
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.