

Product Datasheet

TRAIL-R1 Recombinant Protein (orb1230152)

Description	TRAIL-R1 Recombinant Protein
Target	TNFRSF10A
Preservatives	Contains PBS. Reconstitute with 50 µl sterile water.
Form/Appearance	Lyophilized
Concentration	1mg/ml after reconstitution.
Storage	Stable for at least 6 months after receipt when stored at -20°C.
Tag	Fc Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	The extracellular domain of human TRAIL-R1 (aa 24-239) is fused to the Fc portion of human IgG1.
Purity	> 95% (SDS-PAGE). Endotoxin level is less than 0.01EU/ µg purified protein (LAL test; Lonza).
Source	HEK293 cells
Uniprot ID	Q62009
NCBI	O00220
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.