

Product Datasheet

RANKL Recombinant Protein (orb1230157)

Description	RANKL Recombinant Protein
Target	Tnfsf11
Preservatives	Contains PBS. Reconstitute with 100 μ l sterile water.
Form/Appearance	Lyophilized
Concentration	0.1mg/ml after reconstitution.
Storage	Stable for at least 6 months after receipt when stored at -20°C.
Tag	Fc Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	Mouse RANKL (aa 157-316) is fused at the N-terminus to the Fc portion of human IgG1.
Purity	> 95% (SDS-PAGE). Endotoxin level is less than 0.01EU/ μ g purified protein (LAL test; Lonza).
Source	HEK293 cells
Uniprot ID	O14763
NCBI	O35235
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.