

Product Datasheet

TNF-alpha multimeric Recombinant Protein (orb1230176)

Description	TNF-alpha multimeric Recombinant Protein
Target	TNF
Preservatives	Contains PBS. Reconstitute with 100 µl sterile water.
Form/Appearance	Lyophilized
Concentration	0.1mg/ml after reconstitution.
Storage	Stable for at least 6 months after receipt when stored at -20°C.
Tag	DDDDK Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	The extracellular domain of human TNF-alpha (aa 85-233) is fused at the N-terminus to human ACRP30headless (aa 16-108) and a DDDDK-tag.
Purity	> 95% (SDS-PAGE). Endotoxin level is less than 0.02EU/ µg purified protein (LAL test; Lonza).
Source	HEK293 cells
Uniprot ID	P01584
NCBI	P01375
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.