

Product Datasheet

BAFF Recombinant Protein (orb1230179)

Description	BAFF Recombinant Protein
Target	TNFSF13B
Preservatives	Contains PBS. Reconstitute with 100 µl sterile water.
Form/Appearance	Lyophilized
Concentration	0.1mg/ml after reconstitution.
Storage	Stable for at least 6 months after receipt when stored at -20°C.
Tag	DDDDK Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	Human BAFF (aa 134-285) is fused at the N-terminus to a DDDDK-tag. Natural cleaved form.
Purity	> 95% (SDS-PAGE). Endotoxin level is less than 0.01EU/ µg purified protein (LAL test; Lonza).
Source	E. coli
Uniprot ID	Q9Y6Q6
NCBI	Q9Y275
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.