

Product Datasheet

LIGHT, Soluble Recombinant Protein (orb1230183)

Description	LIGHT, Soluble Recombinant Protein
Target	TNFSF14
Preservatives	Contains PBS. Reconstitute with 100 µl sterile water.
Form/Appearance	Lyophilized
Concentration	0.1mg/ml after reconstitution.
Storage	Stable for at least 6 months after receipt when stored at -20°C.
Tag	DDDDK Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	Human LIGHT (aa 89-240) is fused at the N-terminus to a DDDDK-tag.
Purity	> 95% (SDS-PAGE). Endotoxin level is less than 0.01EU/ µg purified protein (LAL test; Lonza).
Source	CHO cells
Uniprot ID	P29965
NCBI	O43557
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.