
Product Datasheet

TNF-alpha, Soluble Recombinant Protein (orb1230185)

Description

TNF-alpha, Soluble Recombinant

Target	TNF
Preservatives	Contains PBS. Reconstitute with 50 µl sterile water. PBS containing at least 0.1% rAlbumin should be used for further dilutions.
Form/Appearance	Lyophilized
Concentration	10 ug size: 0.1 mg/ml after reconstitution. 50 ug size: 1 mg/ml after reconstitution.
Storage	Stable for at least 1 year after receipt when stored at -20°C.
Tag	DDDDK Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	A DDDDK-tag is fused at the N-terminus of human TNF-alpha (aa 85-223).
Purity	> 95% (SDS-PAGE). Endotoxin level is less than 0.01EU/ µg purified protein (LAL test; Lonza).
Source	E. coli
Uniprot ID	O14788
NCBI	P01375
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.