

Product Datasheet

LSM5 antibody (orb12302)

Description	Goat polyclonal to LSM5.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Target	LSM5
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize isoform a NP_036454.1 only.
Protein Sequence	KCIGSRIHIVMKSDK
Clonality	Polyclonal
MW	9.94; 6.85
Entrez	23658
Dilution Range	ELISA: 1:8000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.