

Product Datasheet

DLL4 Recombinant Protein (orb1230243)

Description	DLL4 Recombinant Protein
Target	DLL4
Preservatives	Contains PBS + 0.5 % Trehalose. Reconstitute with 50 µl sterile water.
Form/Appearance	Lyophilized
Storage	Stable for at least 6 months after receipt when stored at -20°C.
Tag	Fc Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	Signal peptide and extracellular domain of human DLL4 (aa 1-529) are fused at the C- terminus to the Fc portion of human IgG1.
Purity	> 95% (SDS-PAGE). Endotoxin level is less than 0.01EU/ µg purified protein (LAL test; Lonza).
Source	HEK293 cells
Uniprot ID	P55083
NCBI	Q9NR61
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.