

Product Datasheet

IL-33 Recombinant Protein (orb1230270)

Description	IL-33 Recombinant Protein
Target	IL33
Preservatives	0.2um-filtered solution in 55mM TRIS-Cl, pH 8.2, containing 150mM NaCl.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Working aliquots are stable for up to 3 months when stored at -20°C.
Tag	His Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	Human IL-33 (aa 112-270) is fused at the C-terminus to a His-tag.
Purity	> 90% (SDS-PAGE). Endotoxin level is less than 1EU/ µg purified protein (LAL test; Lonza).
Source	E. coli
Uniprot ID	Q8NI99
NCBI	Q2YEJ5
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.