

## Product Datasheet

### Nampt Recombinant Protein (orb1230277)

<b>Description</b>	Nampt Recombinant Protein
<b>Target</b>	NAMPT
<b>Preservatives</b>	Contains PBS. Reconstitute with 100 µl sterile water.
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Working aliquots are stable for up to 3 months when stored at -20°C.
<b>Tag</b>	DDDDK Tag
<b>Note</b>	For research use only
<b>Application notes</b>	This recombinant proteins is for research use only.
<b>Protein Sequence</b>	Extracellular human Nampt (visfatin/PBEF) (aa 1-491) is fused at the N-terminus to a DDDDK-tag.
<b>Purity</b>	> 95% (SDS-PAGE). Endotoxin level is less than 0.01EU/ µg purified protein (LAL test; Lonza).
<b>Source</b>	HEK293 cells
<b>Uniprot ID</b>	<a href="#">Q7TS55</a>
<b>NCBI</b>	<a href="#">P43490</a>
<b>Dilution Range</b>	This recombinant proteins is for research use only.
<b>Expiration Date</b>	6 months from date of receipt.