

Product Datasheet

Nampt Recombinant Protein (orb1230280)

Description	Nampt Recombinant Protein
Target	NAMPT
Preservatives	Salt-free. Reconstitute in distilled water.
Form/Appearance	Lyophilized
Storage	Working aliquots are stable for up to 3 months when stored at -20°C.
Tag	His Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	Intracellular human Nampt (visfatin/PBEF) (aa 1-491) is fused at the C-terminus to a His-tag.
Purity	> 90% (SDS-PAGE). Endotoxin level is less than 1EU/ µg purified protein (LAL test; Lonza).
Source	E. coli
Uniprot ID	Q8R1L8
NCBI	P43490
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.