

Product Datasheet

IL-37 Recombinant Protein (Animal Free) (orb1230377)

Description	IL-37 Recombinant Protein (Animal Free)
Target	IL37
Form/Appearance	Lyophilized
Storage	The recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted protein is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing the recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	MVHTSPKVKN LNPCKFSIHD QDHKVLVLDS GNLIAVPDKN YIRPEIFFAL ASSLSSASAE KGSPILLGVS KGEFCLYCDK DKGQSHPSLQ LKKEKLMKLA AQKESARRPF IFYRAQVGSW NMLESAAHPG WFICTSCNCN EPVGVTDKFE NRKHIEFSFQ PVCKAEMSPS EVSD
Purity	≥ 95% by SDS-PAGE gel and HPLC analyses.
Source	E. coli
Uniprot ID	Q9NZH6
Expiration Date	6 months from date of receipt.