

Product Datasheet

IL-17D Recombinant Protein (Animal Free) (orb1230380)

Description	IL-17D Recombinant Protein (Animal Free)
Target	IL17D
Form/Appearance	Lyophilized
Storage	The recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted protein is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing the recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	APRAGRRPAR PRGCADRPEE LLEQLYGRLA AGVLSAFHHT LQLGPREQAR NASCPAGGRP ADRRFRPPTN LRSVSPWAYR ISYDPARYPR YLPEAYCLCR GCLTGLFGEE DVRFRSAPVY MPTVVLRRTP ACAGGRSVYT EAYVTIPVGC TCVPEPEKDA DSINSSIDKQ GAKLLLGPND APAGP
Purity	≥ 98% by SDS-PAGE gel and HPLC analyses.
Source	E. coli
Uniprot ID	Q8TAD2
Expiration Date	6 months from date of receipt.